

## VARIABLE LOUPES

### **LPM-P35V-4407**

Zoom magnification **2.5x to 3.5x**  
Variable working distance 440mm-250mm  
Variable depth of field 130mm-50mm  
Variable field of view 100mm-40mm



Weight 2.5oz



## 3.5X PREMIUM LOUPES



Weight 2.5oz

SHORT: Close-up, bench work

### **LPM-P35F-4007 (3.5x)**

-Working distance 365 +/- 15mm  
-Depth of field 30mm  
-Field of view 44mm

LONG: Stand-up surgeries

### **LPM-P35F-4607 (3.5x)**

-Working distance 460 +/- 20mm  
-Depth of field 40mm  
-Field of view 45mm

## 2.5X PREMIUM LOUPES

SHORT: Close-up, bench work

### **LPM-P25F-3007 (2.5x)**

-Working distance 340 +/- 25mm  
-Depth of field 50mm  
-Field of view 66mm

LONG: Stand-up surgeries

### **LPM-P25F-4607 (2.5x)**

-Working distance 460 +/- 30mm  
-Depth of field 80mm  
-Field of view 90mm



Weight 2.6oz



### **LPM-C25F-3807 (2.5x)**

-Working distance 380 +/- 30mm  
-Depth of field 60mm  
-Field of view 80mm

### **LPM-C35F-3807 (3.5x)**

-Working distance 380 +/- 20mm  
-Depth of field 40mm  
-Field of view 50mm

Works on most prescription glasses

## CLIP-ON LOUPES



Weight 2.3 oz

